

Instruction Manual

Folding Aluminium Walker with Seat



FWALSEATMSIL

SectorCare

<https://sectorcare.com.au>



Product Overview



1 Hand Grip

2 Push-Button
Folding Mechanism

3 Main Frame

4 Seat

5 Push-Button Folding
Mechanism

6 Height Adjustable
Leg

Product Installation

Parts List:



Main Frame*1



Seat*1



Hex Screw*2



Hex Nut*2

Take all items out of the box.



Unfold the main frame as shown in the figure.



Place the seat board flat on the main frame, press the buckles to lock and secure, as shown in the figure.



Insert hex screws and hex nuts into the buckle positions on both sides of the seat board respectively, then tighten them for fixation, as shown in the figure.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.