

Instruction Manual

Folding Aluminium Walker with Wheels and Brakes



FWALWBMSIL

SectorCare

Product Overview



1 Hand Grip

2 Frame Folding Button

3 Height Adjustable
Leg

4 Front Wheels

5 Anti-slip Rubber Tips

Product Installation

1. Take out the walker and all of the accessories from the carton.
2. Swing both side frames outward from the front cross brace and unfold frame completely until an audible "click" is heard.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

SectorCare



<https://sectorcare.com.au>