

Instruction Manual

Folding Aluminium Walker with Wheels



FWALWMSIL

SectorCare

<https://sectorcare.com.au>



Product Overview



1 Hand Grip

2 Push-Button
Folding Mechanism

3 Main Frame

4 Height Adjustable Leg

5 Anti-slip Rubber Tips

6 Front Wheel

Product Installation

Parts List:



Main Frame*1



Front Wheel*2




Foot Tube*2

Open the carton and take out the product, then unfold the main frame as shown in the figure.



Align the 2 foot tubes and 2 wheeled foot tubes with the corresponding connectors respectively, insert and secure them as shown in the figure.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.