

Instruction Manual

Bath Chair with Handle and Back 360 Swivel



SCAL36001SSIL

SectorCare

Product Overview



- 01 Seat
- 02 Armrest
- 03 Height Adjustable Leg
- 04 Non-slip Rubber Feet

Assembly Instructions



Seat x1



Base x1



Instruction Manual x1

Product Installation

1. Remove the product from the packaging.



2. Place the seat pad on top of the base frame, align with the central cylindrical interface, insert and press down firmly to secure, as illustrated:



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.

SectorCare



<https://sectorcare.com.au>