

Instruction Manual

Bath Chair with Handle and Back 360 Swivel No Legs



SCST36001WHI

SectorCare

<https://sectorcare.com.au>



Product Installation



01 Seat Board

02 Armrest Frame

03 360° Rotatable Handle

Parts List:



Seat Board*1



Base Support Frame*1

1. Take out all products and accessories after opening the package.



2. Place the seat flat on top of the base, align it with the central cylindrical connector of the base, insert and press down to secure, as shown in the figure.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.