

Instruction Manual

Toilet Aid with Armrest



THRAL01SSIL

SectorCare

<https://sectorcare.com.au>



Product Installation



01 Armrest

02 Toilet Seat

03 Adjustable Knob

Parts List:



Armrest Grip*2



Commode Seat*1

1. Take out all products and accessories after opening the package.



2. Align the armrest with the grooves on both side fittings of the toilet, lay flat and press down firmly. Place nuts into the matching positions on the white plastic section of the armrest and secure tightly, as shown in the figure.



3. Insert screws into the connectors on both sides of the toilet seat, then tighten to secure, as shown in the figure.



LEGAL STATEMENT

This manual and its contents are protected under the Copyright Act 1968 (Cth). No part of this publication may be reproduced, distributed, modified, stored in a retrieval system, or transmitted in any form or by any means without the prior written permission of FRG Pty Ltd. Unauthorized use is strictly prohibited.

© 2026 FRG Pty Ltd. All rights reserved.